

## **Frozen Tissue Xtractor:**

**Preserves sample quality and streamlines sample processing.**

- Simple and easy design for extracting out small aliquots of tissue.
- Provides smaller labs workflow options.
- Supports highest sample quality and integrity without the **THAW** cycle.
- Preserves frozen tissue without contamination allowing for increase sample utilization.



### **Kits comes with:**

**Handheld Tissue Xtractor, 6 bags of coring probes (probes), metal tissue tray holder, probe chilling nest and a dry ice bucket.**

**Optional items: Cryostat adapter with plastic tissue trays.**

## **Benefits of The Tissue Xtractor**

**Eliminates the **thaw** cycle preserving labile molecules**

**Various tissue types and sizes can be cored**

**Targets specific tissue sections.**

**Reduces operating costs, increasing productivity**

**Avoids manual slicing of tissue**

**Extracts out small volume aliquots for downstream analysis**



# FROZEN ALIQUOTTING SOLUTIONS